

CHAPTER 5

NIRT MONITORING PROGRAM

Section Number	Section	Page Number
5.1	GENERAL INFORMATION .....	5-1
5.2	REVIEW OF SRS WORKSHEETS .....	5-2
5.3	CHECK SAMPLES .....	5-2
5.4	MONITORING FILE SAMPLES .....	5-3
5.5	WHEAT/BARLEY/SOYBEAN QUALITY CONTROL CHARTS .....	5-5
5.6	COLLABORATIVE CHECK SAMPLES AND SPECIAL STUDIES .....	5-9
5.7	INTERMARKET SAMPLE EXCHANGE .....	5-10
	Figure 1 – Average Difference Between Field Location and TSD .....	5-11
	Figure 2 – Range Difference Between Field Location and TSD .....	5-12



CHAPTER 5

NIRT MONITORING PROGRAM

**5.1 GENERAL INFORMATION**

- a. The NIRT monitoring program is designed to monitor the accuracy of official wheat, barley, soybean, and corn testing services. Wheat wet gluten content is determined directly from the NIRT protein, and no additional monitoring is required. TSD will review calibration performance using the existing HRW and HRS protein monitoring sample submissions.
- b. Accuracy is evaluated by comparing field office and specified service point results with TSD results. Field offices with NIRT instruments may conduct supplemental file sample monitoring of specified service points within the circuit in addition to TSD monitoring to evaluate the performance of instruments and technicians.
- c. Monitoring information is used by specified service points to evaluate the performance of their testing program. Field offices use the monitoring information to evaluate their own accuracy as well as the accuracy of the testing service program within their circuit. TSD uses the monitoring information to evaluate the accuracy of individual field locations as well as the entire program.
- d. TSD monitoring will identify service points generating questionable results. Field offices and specified service points must initiate corrective action and follow up whenever wheat, barley, soybean, or corn testing problems are detected. Corrective action and follow-up includes investigating, identifying, correcting, and documenting the cause of accuracy problems.
- e. Several methods are utilized to monitor particular elements in the testing process. These methods include review of SRS worksheets, check samples, monitoring file samples (quality control charts), intermarket sample exchanges, and special studies. Procedures regarding these monitoring methods are discussed in the following sections.

## 5.2 REVIEW OF SRS WORKSHEETS

TSD, field offices, and specified service points record the results of daily and weekly standard reference sample checks on SRS worksheets (examples are shown in Chapter 3).

- a. Frequency. Locations are no longer required to submit a copy of all SRS information recorded for each weekly period to the monitoring field office. TSD, or the monitoring field office manager may request SRS information if a problem is suspected (unusual high or low results, board appeal, foreign compliant, etc.) NIRT coordinators and specified service points must review these worksheets for accuracy and completeness.
- b. Evaluation of Results. NIRT coordinators must evaluate the SRS worksheets to determine if: (1) bias adjustments were completed when necessary, (2) instrument accuracy was maintained, or (3) bias adjustments were required frequently (which may indicate the need for operator training, instrument repair, or SRS replacement).

## 5.3 CHECK SAMPLES

- a. Corn. Check samples are used to test the capability of the national corn protein, oil, and starch testing program. Corn check samples are issued once per year (prior to harvest) by TSD. Specified service points and laboratories that provide corn protein, oil, and starch testing services for FGIS must participate in the check sample program.
  - (1) Check samples originate from TSD and will include specific instructions for testing the samples. Upon receipt, inspection points must complete the analysis of these samples within 10 working days. Report check sample results to TSD on the forms provided. Whenever possible, the results should be faxed to TSD to expedite the preparation of the check sample report.
  - (2) The field office will be notified immediately if the test results indicate a problem at a test location. Otherwise, TSD will not release any information on the content of the check samples until the results from all laboratories have been received. TSD will analyze the data and issue a report which will be distributed to all participating official agencies and supervising field offices.
- b. Soybeans. The soybean protein and oil check sample program has been discontinued. Soybean protein and oil accuracy are now evaluated by the file sample monitoring program and other methods.

#### 5.4 MONITORING FILE SAMPLES

Wheat, barley, and soybean testing accuracy are evaluated through the file sample monitoring system. Official inspection points select and forward samples to TSD. Accuracy is determined by comparing original results with the average of the TSD master instruments' results. (If a field office has a local, supplemental monitoring program, official agencies shall also forward samples to the field office location.)

a. Selecting Samples.

- (1) Select monitor samples representing the range of constituents observed during the week. Include a low and a high protein sample each week with three intermediate samples. When applicable, include a high oil and low oil sample.
- (2) Avoid selecting all samples tested from the same day. When less than five samples of a wheat class or barley subclass are tested during a week, select all samples tested for that class or subclass. Do not make up monitoring samples to fulfill the minimum number of samples. Do not select mixed wheat samples.
- (3) Each dockage-free wheat or barley sample submitted for monitoring must be at least 600 grams unless special arrangements have been made with TSD. Each soybean sample should be at least 800 grams.

b. Select five samples of wheat per week. Do not submit five samples for each wheat class tested; submit five samples total. The samples should all be the same class. Locations routinely certifying NIRT protein on more than one wheat class should, where practical, rotate weekly submissions among those classes.

c. Select five samples per barley subclass tested, per week.

d. Select five samples of soybeans tested, per week.

e. Field offices and/or TSD may request additional samples for monitoring and/or special study purposes.

f. Mailing Instructions. Seal the samples in individual 6-mil plastic bags. Mark the sample number and wheat class or barley subclass on each bag using an indelible marker.

Place the samples and sample information (scan sheet) in a canvas mailing bag and indicate "NIRT Monitoring Sample" on the reverse side of the mail tag (use an orange color mailing tag). This will assist in separating protein monitoring samples from other samples received in the mail. According to postal requirements, Federal employees must use metered mailing tags. All nonfederal employees (i.e., delegated and designated agencies) must use business reply mailing tags.

- g. Monitoring Results. Field office personnel must retrieve the monitoring data from the Quality Assurance/Quality Control (QAQC) homepage—<http://nqdbweb.gipsa.usda.gov:8010/>, review the information, and immediately initiate follow-up action when accuracy deficiencies are indicated by monitoring results. In addition, field office managers must forward TSD monitoring results to the specified service points for their review and follow-up action within 2 working days of retrieval.
- h. Evaluating Monitoring Results. Field offices and specified service points will evaluate completed quality control charts to determine if any action limit (tolerance limits, absolute limits, and/or run limit) violations occurred.
  - (1) Action limit violations occurring on the average difference chart generally indicate a bias-related problem. Action limit violations occurring on the range difference chart generally indicate inconsistency due to fluctuating laboratory conditions, failure to follow procedures consistently, instrument problems, or improper instrument slope or bias adjustment. Violations on the range difference chart are actually more serious than those on the average difference chart because if the results are inconsistent, the average differences are not meaningful.
  - (2) Monitoring field offices and agencies must initiate corrective action when quality control chart rule violations occur. Field office managers must document any action taken to resolve the differences. This documentation includes action taken to identify the cause and extent of the problem and steps to resolve the problem and/or reasons why no further action is necessary. Documentation may be placed directly on the control chart indicating action limit violations.

## 5.5 WHEAT/BARLEY/SOYBEAN QUALITY CONTROL CHARTS

a. General. A Quality Control Chart (QCC) is a visual display of monitoring data. A QCC effectively displays extreme variations, shifts, and trends. Also, a QCC illustrates the difference between results while statistically defining expected variability using control limits. These limits are established based on the normal expected variation of results.

b. Control Charts. The protein/oil monitoring program utilizes Average Difference and Range Difference quality control charts. The average difference chart illustrates the difference between a specified service point's average for five weekly monitoring samples and the TSD's average for the same samples. The range difference chart plots the range of individual sample differences for the corresponding weekly monitoring sample set.

(1) Average Difference Chart. The average difference chart includes a zero or Center Line (CL), upper and lower Tolerance Limits (TL), and upper and lower Absolute Limit lines (AL).

The center line is the control chart reference point. Points plotted above the center line indicate a positive difference when compared to the monitoring result. Points plotted below the center line indicate a negative difference when compared to the monitoring result.

- (a) For wheat, the absolute limit lines are set at  $\pm 0.20$  percent protein from the center line. The tolerance limit lines are set at  $\pm 0.15$  percent protein from the center line. The run limit is set at  $\pm 0.10$  percent protein from the center line.
- (b) For barley, the absolute limit lines are set at  $\pm 0.25$  percent protein from the center line. The tolerance limit lines are set at  $\pm 0.20$  percent protein from the center line. The run limit is set at  $\pm 0.10$  percent protein from the center line.
- (c) For soybean protein, the absolute limit lines are set at  $\pm 0.25$  percent protein from the center line. The tolerance limit lines are set at  $\pm 0.20$  percent protein from the center line. The run limit is set at  $\pm 0.15$  percent protein from the center line.

- (d) For soybean oil, the absolute limit lines are set at  $\pm 0.20$  percent oil from the center line. The tolerance limit lines are set at  $\pm 0.15$  percent oil from the center line. The run limit is set at  $\pm 0.10$  percent oil from the center line.

These tolerances are determined statistically based on the systems actual performance and may be revised from time to time. A violation of any of the established tolerances means that there is less than one chance in one hundred that the observed error level occurred due to random chance. Therefore, it is very likely that a correctable problem exists. Average and Range Difference control charts are shown in Figures 1 and 2.

- (2) Range Difference Chart. The range difference chart indicates how much difference is observed within a set of monitoring samples.
  - (a) For wheat, the absolute limit line is set at 0.60 percent protein. The tolerance limit line is set at 0.40 percent protein.
  - (b) For barley, the absolute limit line is set at 0.70 percent protein. The tolerance limit line is set at 0.50 percent protein.
  - (c) For soybeans, the absolute limit line is set at 0.80 percent protein. The tolerance limit line is set at 0.40 percent protein.
  - (d) For soybeans, the absolute limit line is set at 0.60 percent oil. The tolerance limit line is set at 0.45 percent oil.

Table 1. Action Limits

Violation If	Average Difference Chart			Range Difference Chart	
	Avg. diff. exceeds limit on any run	Avg. diff. exceeds limit on 2 consecutive runs	(a) 4 of 5 consecutive runs are all either + or -, (b) all 4 runs exceed limit	Range of differences exceeds limit	Range of differences exceed limit on 2 consecutive runs
	Absolute Limit	Tolerance Limit	Run Limit	Absolute Limit	Tolerance Limit
Wheat Protein	0.20	0.15	0.10	0.60	0.40
Barley Protein	0.25	0.20	0.10	0.70	0.50
Soybean Protein	0.25	0.20	0.15	0.80	0.60
Soybean Oil	0.20	0.15	0.10	0.60	0.45

c. Plotting Control Charts. Control charts are generated as follows.

- (1) Sequential Plotting. Control charts depict how a process or system is varying over time. In order for the results to be meaningful, data is plotted in the sequence in which the original results were obtained. That is, when plotting the Average/Range chart, samples originally tested during week 3 are plotted after samples originally tested during week 2, but before samples originally tested during week 4 (see **Figures 1 and 2**).
- (2) Average Difference. The difference between the original average result and the monitoring average results is obtained by subtracting the monitoring office result from the original protein result (see **example**).

- (3) Range Difference. The difference between the original result and the monitoring result is determined for each individual sample monitored within a sample set. The range difference value that is plotted is the difference between the smallest and the largest difference observed in one set. If these two differences are of the opposite sign, add the magnitudes of the two numbers. The range difference is calculated and plotted on the range difference chart.
- (4) Plotted Results. The difference between the monitoring average result and the original result obtained from step 2 is plotted on the average difference chart. The positive differences (original results higher in protein than the monitoring results) are plotted above the center line. Negative differences (original results lower in protein than the monitoring results) are plotted below the center line. Results having no difference (original result is identical to monitoring result) are plotted on the center line.

Example of Plotting Control Chart				
Week Ending Date	Sample No.	Original Protein	TSD Protein	ORIG-TSD Difference
6/28/99	1	13.15	13.07	+0.08
6/28/99	2	10.18	10.54	-0.36
6/28/99	3	14.47	14.79	-0.32
6/28/99	4	12.56	12.23	+0.33
6/28/99	5	18.57	18.77	-0.20
Average		13.79	13.88	-0.09

In the example, the average difference between the service point and TSD is -0.09.

The range difference is calculated by determining the largest spread between the individual sample differences. In this example, the largest difference is +0.33 for sample No. 4 and -0.36 for sample No. 2. The spread or range between these two numbers is 0.69. Therefore, 0.69 is the range difference. Figures 1 and 2 illustrate average and range difference charts.

- d. Reviewing Control Charts—Action Limits. Protein testing problems are indicated on the control charts either by a large difference between the average protein results or a consistent pattern of smaller differences on a series of average results. Three action limits are used for rapid identification of protein testing problems through interpretation of the control chart.
- (1) Absolute Limit (AL). This action limit is intended to identify excessive differences between results and indicates a potential protein testing problem. An absolute limit violation occurs if any plotted value is equal to or greater than the absolute limit lines.
  - (2) Tolerance Limit (TL). This action limit controls the number of consecutive data sets with a large difference between original and monitoring results but not so large as to exceed the absolute limit. An average difference limit violation occurs if two consecutive data sets are either equal to or above the upper tolerance limit line or both equal to or below the lower tolerance limit line. A range difference tolerance limit violation occurs if two consecutive data sets are equal to or greater than the tolerance limit line.
  - (3) Run Limit. This action limit controls the number of consecutive comparisons which are all above or all below the average difference chart center line (CL). A violation occurs if four out of five consecutive results are either all above or all below the CL, and the average difference from the center line for these four results exceed the applicable limit from Table 1. Run limits do not apply to the range difference chart.

## 5.6 COLLABORATIVE CHECK SAMPLES AND SPECIAL STUDIES

- a. Collaborative Check Samples. Collaborative check samples may be initiated by TSD for cross-checking other data. TSD will select and send enough of each sample selected so that specified service points can retain a portion for rechecking purposes. Upon receipt, participants must complete the analysis and report the data within 5 working days.

Use the forms provided with the samples to report analysis results. Retain a copy of the completed form and return the original immediately to TSD. Results from all locations included in the collaborative study will be compiled by TSD and reported to participating field offices and specified service points.

- b. Special Studies. TSD may initiate and conduct special studies. These studies are designed for a specific purpose (i.e., resolving differences either within or between markets, evaluating calibration performance, or updating the calibrations). When special studies are initiated by TSD, it is required that all participants (as designated by TSD) respond with utmost priority, as these are normally of an urgent nature and an expedient resolution of the problem is essential. Because TSD will not be routinely monitoring corn file samples, the field office will be required to provide TSD with samples for monitoring and updating the calibration. Periodically throughout the year, TSD will contact selected field offices in corn markets and request that 10 randomly selected corn samples be provided to monitor calibration performance.


During harvest, TSD will request those field locations providing corn NIRT services to collect samples within specified ranges of protein, oil, and starch for use in updating the calibrations. Prior to harvest, TSD will notify the field office of the types and numbers of samples needed and instructions for shipping them to TSD.

## **5.7 INTERMARKET SAMPLE EXCHANGE**

An intermarket sample exchange helps isolate protein and/or oil differences between inspection points. Protein and oil testing laboratories will determine protein and oil results on separate portions obtained from the same sample. Protein and oil results are then compared to determine if any significant differences exist.

There are no restrictions as to which offices may exchange samples. Specified service points are encouraged to exchange samples with other specified service points and field offices for the purpose of resolving intermarket inspection differences. A copy of the results of the exchange must be provided to the field office and/or TSD for review if they were not participants in the exchange.

Figure 1



AL = Absolute Limit

CL = Center Line

\* = Average of 5 samples


TL = Tolerance Limit

1, 2, 3, or 4 = Number  
samples averaged

Figure 2

### Range Difference Between Field Location and TSD — HRW

SERVICE POINT FULL NAME= CULVERTON CLASS=HRW



AL = Absolute Limit

TL = Tolerance Limit

### **EEO STATEMENT**

The U.S. Department of Agriculture (USDA) prohibits discrimination in its programs on the basis of race, color, national origin, sex, religion, age, disability, political beliefs, and marital or familial status. (Not all prohibited bases apply to all programs.) Persons with disabilities who require alternate means for communication of program information (Braille, large print, audiotape, etc.) should contact USDA's TARGET Center at (202) 720-2600 (voice and TDD). To file a complaint, write to the USDA, Office of Civil Rights, Room 326-W, 1400 Independence Avenue, SW, Washington, DC 20250-9410, or call (202) 720-5964 (voice and TDD). USDA is an equal employment opportunity employer.